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### TURKEY COUNTRY SPECIFIC INFORMATION:

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### BRAZIL COUNTRY SPECIFIC INFORMATION:

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### INDICATIONS FOR USE (NON-EU REGIONS)

The Irregular Rhythm Notification Feature is a software-only mobile medical application that is intended to be used with the Apple Watch. The feature analyses pulse rate data to identify episodes of irregular heart rhythms suggestive of atrial fibrillation (AF) and provides a notification to the user. The feature is intended for over-the-counter (OTC) use. It is not intended to provide a notification on every episode of irregular rhythm suggestive of AF and the absence of a notification

is not intended to indicate no disease process is present; rather the feature is intended to opportunistically surface a notification of possible AF when sufficient data are available for analysis. These data are only captured when the user is still. Along with the user's risk factors, the feature can be used to supplement the decision for AF screening. The feature is not intended to replace traditional methods of diagnosis or treatment.

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### **Set-Up/On-boarding**

- For region availability and device compatibility for the IRNF, please visit <https://support.apple.com/HT208931>
- Update Apple Watch and iPhone to latest OS.
- Open the Health app on your iPhone and select "Browse".
- Navigate to "Heart", then select "Irregular Rhythm Notifications".
- Follow the on-screen instructions.
- You may exit on-boarding at any time by tapping "Cancel".

### **Receiving a Notification**

- Once the feature is turned on, you will receive a notification if the feature identified a heart rhythm suggestive of AF and confirmed it on multiple readings.
- If you have not been diagnosed with AF by a physician, you should discuss the notification with your GP.

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New data cannot be collected once your Apple Watch's storage is full. You should free up space by deleting unwanted apps, music or podcasts. You can check your storage usage by navigating to the Apple Watch app on your iPhone, tapping "My Watch", tapping "General" and then tapping "Usage".

## **SAFETY AND PERFORMANCE**

In a study of 226 participants aged 22 years or older who had received an AF notification while wearing Apple Watch and subsequently wore an electrocardiogram (ECG) patch for approximately 1 week, 41.6% (94/226) had AF detected by ECG patch. During concurrent wear of Apple Watch and an ECG patch, 57/226 participants received an AF notification. Of those, 78.9% (45/57) showed concordant AF on the ECG patch and 98.2% (56/57) showed AF and other clinically relevant arrhythmias. A total of 370 irregular rhythm notifications with readable ECG patch data was received by the 57 participants. Of those 370 notifications, 322 (87.0%) were assessed to be AF, 47 (12.7%) were arrhythmias other than AF, and 1 (0.3%) was sinus rhythm. These results demonstrate that, while in the majority of cases the notification will accurately represent the presence of AF, in some instances, a notification may indicate the presence of an arrhythmia other than AF. No serious device adverse effects were observed.

## CAUTIONS

**The Irregular Rhythm Notification Feature cannot detect heart attacks. If you ever experience chest pain, pressure, tightness or what you think is a heart attack, call emergency services.**

**The Irregular Rhythm Notification Feature is not constantly looking for AF and should not be relied on as a continuous monitor. This means the feature cannot detect all instances of AF and people with AF may not get a notification.**

Apple Watch may be unable to collect data when Apple Watch is in close vicinity to strong electromagnetic fields (e.g. electromagnetic anti-theft systems, metal detectors).

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DO NOT wear your Apple Watch during a medical procedure (e.g., magnetic resonance imaging, diathermy, lithotripsy, cautery and external defibrillation procedures).

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Not intended for use by individuals under age 22.

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Apple does not guarantee that you are not experiencing an arrhythmia or other health conditions even in the absence of an irregular rhythm notification. You should notify your physician if you experience any changes to your health.

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**SECURITY:** Apple recommends that you add a passcode (personal identification number [PIN]), Face ID or Touch ID (fingerprint) to your iPhone and a passcode (personal identification number [PIN]) to your Apple Watch to add a layer of security. It is important to secure the iPhone since you will be storing personal health information.

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---

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ENGLISH (UK)

## Irregular Rhythm Notification Feature

### Instructions for Use

Apple Inc.  
One Apple Park Way  
Cupertino, CA 95014  
[www.apple.com](http://www.apple.com)

Apple Distribution International  
Hollyhill Industrial Estate,  
Hollyhill, Cork,  
Ireland  
Contact: [medicalcompliance@group.apple.com](mailto:medicalcompliance@group.apple.com)

### TURKEY COUNTRY-SPECIFIC INFORMATION:

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Product Tracking Number for Irregular Rhythm Notification Feature: 8682450178007

### BRAZIL COUNTRY-SPECIFIC INFORMATION:

ANVISA N<sup>o</sup>: 80117580881

TECHNICAL DIRECTOR: Luiz Levy Cruz Martins - CRF-SP: 42415

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### INDICATIONS FOR USE (NON-EU REGIONS)

The Irregular Rhythm Notification Feature is a software-only mobile medical application that is intended to be used with Apple Watch. The feature analyses pulse rate data to identify episodes of irregular heart rhythms suggestive of atrial fibrillation (AF) and provides a notification to the user. The feature is intended for over-the-counter (OTC) use. It is not intended to provide a notification of every episode of irregular rhythm suggestive of AF, and the absence of a notification is not

intended to indicate no disease process is present; rather, the feature is intended to opportunistically display a notification of possible AF when sufficient data is available for analysis. This data is only captured when the user is still. Along with the user's risk factors, the feature can be used to supplement the decision for AF screening. The feature is not intended to replace traditional methods of diagnosis or treatment.

The feature has not been tested for and is not intended for use in people under 22 years of age. It is also not intended for use in individuals previously diagnosed with AF.

## **INTENDED PURPOSE (EU REGION)**

### Intended Use

The Irregular Rhythm Notification Feature (IRNF) is intended to pre-screen and notify the user of the presence of irregular rhythms suggestive of atrial fibrillation (AF). The feature can be used to supplement a clinician's decision to screen for possible AF. The feature is intended for over-the-counter (OTC) use.

The feature has not been tested for and is not intended for use in people under 22 years of age. It is also not intended for use in individuals previously diagnosed with AF.

### Indications

The feature is indicated to pre-screen for irregular rhythms suggestive of AF for anyone aged 22 years and over.

## **USING THE IRREGULAR RHYTHM NOTIFICATION FEATURE**

### **Setup/Onboarding**

- For region availability and device compatibility for the IRNF, please visit <https://support.apple.com/HT208931>
- Update Apple Watch and iPhone to latest OS.
- Open the Health app on your iPhone and select Browse.
- Navigate to Heart, then select Irregular Rhythm Notifications.
- Follow the onscreen instructions.
- You may exit onboarding at any time by tapping Cancel.

### **Receiving a Notification**

- Once the feature is turned on, you will receive a notification if the feature identified a heart rhythm suggestive of AF and confirmed it on multiple readings.
- If you have not been diagnosed with AF by a medical practitioner, you should discuss the notification with your doctor.

All data collected and analysed by the Irregular Rhythm Notification Feature is saved to the Health app on your iPhone. If you choose to, you can share that information by exporting your health data in the Health app.

New data cannot be collected once your Apple Watch storage is full. You should free up space by deleting unwanted apps, music or podcasts. You can check your storage usage by navigating to the Apple Watch app on your iPhone, tapping My Watch, tapping General, then tapping Usage.

## **SAFETY AND PERFORMANCE**

In a study of 226 participants aged 22 years or older who had received an AF notification while

wearing Apple Watch and subsequently wore an electrocardiogram (ECG) patch for approximately 1 week, 41.6% (94/226) had AF detected by the ECG patch. During concurrent wear of Apple Watch and an ECG patch, 57/226 participants received an AF notification. Of those, 78.9% (45/57) showed concordant AF on the ECG patch and 98.2 % (56/57) showed AF and other clinically relevant arrhythmias. A total of 370 irregular rhythm notifications with readable ECG patch data were received by the 57 participants. Of those 370 notifications, 322 (87.0%) were assessed to be AF, 47 (12.7%) were arrhythmias other than AF, and 1 (0.3%) was sinus rhythm. These results demonstrate that, while in the majority of cases the notification will accurately represent the presence of AF, in some instances, a notification may indicate the presence of an arrhythmia other than AF. No serious device adverse effects were observed.

## CAUTIONS

**The Irregular Rhythm Notification Feature cannot detect heart attacks. If you ever experience chest pain, pressure, tightness, or what you think is a heart attack, call emergency services.**

**The Irregular Rhythm Notification Feature is not constantly looking for AF and should not be relied on as a continuous monitor. This means the feature cannot detect all instances of AF, and people with AF may not get a notification.**

Apple Watch may be unable to collect data when it is in close proximity to strong electromagnetic fields (e.g. electromagnetic anti-theft systems or metal detectors).

A number of factors can impact the ability of the feature to measure your pulse and detect an irregular rhythm suggestive of AF. These include factors like motion, hand and finger movements, dark tattoos on the wrist and the amount of blood flow to your skin (which can be reduced by cold temperatures).

DO NOT wear your Apple Watch during a medical procedure (e.g. magnetic resonance imaging, diathermy, lithotripsy, cautery and external defibrillation procedures).

DO NOT change your medication without talking to your doctor.

Not intended for use by individuals under the age of 22.

Not intended for use by individuals previously diagnosed with AF.

Notifications made by this feature are potential findings, not a complete diagnosis of cardiac conditions. All notifications should be reviewed by a medical professional for clinical decision-making.

Apple does not guarantee that you are not experiencing an arrhythmia or other health conditions even in the absence of an irregular rhythm notification. You should notify your doctor if you experience any changes to your health.

For best results, make sure your Apple Watch fits snugly on top of your wrist. The heart rate sensor should stay close to your skin.

**SECURITY:** Apple recommends that you add a passcode (personal identification number [PIN]), Face ID or Touch ID (fingerprint) to your iPhone and a passcode (personal identification number [PIN]) to your Apple Watch to add a layer of security. It is important to secure your iPhone since you will be storing personal health information on it.

## EQUIPMENT SYMBOLS

Manufacturer

Read instructions before use

**EC**

**REP**

European Authorised Representative

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